




1

Getting the Most out of this Course

- Ask questions!
- Experiment with hands-on exercises to understand the methods used.
- Explore solutions. Remember that implementations explore one possible solution. You may find a better one!

Course Learning Map



ni.com/training

7

Lesson 1 Moving Beyond Dataflow

TOPICS

- A. Asynchronous Communication
- B. Queues
- C. Event-Driven Programming



ni.com/training

11

A. Asynchronous Communication



| ni.com/training

12

Asynchronous Communication

LabVIEW is a dataflow language.

- Functions depend on other functions for data.
- Dependent functions do not execute until their dependencies have completed execution.
- Wires transfer data between functions.

However, sometimes you need break dataflow and program using asynchronous communication.



| ni.com/training

13

Asynchronous Communication



Asynchronous Communication — Transfer information without using wires

- Asynchronous communication is between the following:
 - Parallel loops
 - UI and Block Diagram
 - VIs
 - Application instances (LabVIEW projects, executables, etc.)
- Information includes the following:
 - Data
 - Notification that something happened



ni.com/training

14

B. Queues

Queues

Queue Operations

Producer/Consumer (Data) Design Pattern



ni.com/training

15

Queues

- Are used for communicating data between parallel loops.
- Store multiple pieces of data (i.e. buffer data).
- Work in a FIFO (first in, first out) manner by default.
- Can hold data of any type.




ni.com/training

16

Drawbacks of Variables

Drawbacks of using variables to communicate between loops:

- It's possible to read duplicate data.
- It's possible to miss data.
- You can create read-modify-write race conditions.



ni.com/training


17

Queue Operations



Use the Queue Operations functions to create and use a queue for communicating data between:


- Different sections of a VI.
- Different VIs.



ni.com/training

18

Producer/Consumer Design Pattern (Data)



ni.com/training

19

C. Event-Driven Programming

Events – Definition

Event-Driven Programming – Definition

Polling Versus Event Structures

Parts of an Event Structure

Configuring the Event Structure

Caveats and Recommendations



| ni.com/training

23

Events





Event — An asynchronous notification that something has occurred

- Events originate from the user interface, external I/O, or other parts of the program.
- Events do things TO event sources.
 - Example: Value change happens TO front panel controls.

24

Event-Driven Programming



| ni.com/training

25

Polling versus Event Structures

Polling

- Method of event-based programming where a loop must continually run code to check if changes have occurred.
- Polling the front panel requires a significant amount of CPU time.
- Polling can fail to detect changes if they occur too quickly.

Event Structures

- Events in Event structures eliminate the need to poll the front panel.
- Benefits of using Event structures:
 - Reduces the CPU requirements of the program.
 - Simplifies the block diagram code.
 - Guarantees that the block diagram can respond to all interactions the user makes.



| ni.com/training

26

Using Event Structures for Event-Driven Programming

An Event structure works like a Case structure with a built-in Wait on Notification function.

Use an Event structure to handle user-interface (static) events such as:

- Pressing a button on the mouse.
- Pressing a key on the keyboard.
- Changing the value of a numeric control.




ni.com/training

27

Parts of an Event Structure


- Event Selector Label—Identifies the event case viewed.
- Timeout—Specifies time in ms to wait for events. Default value is -1 (indefinite).



28

Parts of an Event Structure (continued)

- Event Data Node—Identifies the data LabVIEW provides when the event occurs; similar to the Unbundle By Name function.
- Event Filter Node—Identifies the subset of data available in the Event Data node that the event case can modify.




ni.com/training

29

Using an Event Structure

In general, place Event structures inside While Loops.


- Event structures handle exactly one event per iteration of the While Loop.
- Event structures sleep when no events occur.



ni.com/training

30

Configuring the Event Structure




Use a dialog box to configure each event by right-clicking the Event structure border and selecting **Edit Events Handled by This Case** from the shortcut menu.



ni.com/training

31


Edit Events Dialog Box



ni.com/training

32

Notify and Filter Events



➡ Notify Events (green arrow)

User action has already occurred and LabVIEW has processed the event.

➡ Filter Events (red arrow)

User action has already occurred and LabVIEW has NOT processed the event. Filter events allow you to override default behavior for event.



ni.com/training

33

Caveats and Recommendations

- Avoid using an Event structure outside of a loop.
- Place only one Event structure in a loop.
- Avoid configuring two Event structures for the same event.
- Use a Value Change event to detect value changes.
- Keep event handling code short and quick.
- Place Boolean control terminals inside an event case for latched operations to work properly.



ni.com/training

37

Summary—Quiz

1. Which of the following buffer data?

- a) Queues
- b) Events
- c) Local Variables

38



| ni.com/training

Summary—Quiz Answer

1. Which of the following buffer data?

- a) Queues**
- b) Events
- c) Local Variables

39



| ni.com/training

Summary—Match the Following

1. Obtain Queue



2. Get Queue Status



3. Release Queue



4. Enqueue Element



a. Destroys the queue reference

b. Assigns the data type of the queue

c. Adds an element to the back of a queue

d. Determines the number of elements currently in the queue



ni.com/training

40

Summary—Match the Following Answer

1. Obtain Queue



2. Get Queue Status



3. Release Queue



4. Enqueue Element



a. Destroys the queue reference

b. Assigns the data type of the queue

c. Adds an element to the back of a queue

d. Determines the number of elements currently in the queue



ni.com/training

41

Summary—Quiz

3. Which of the following are valid data types for queues?

- a) String
- b) Numeric
- c) Enum
- d) Array of Booleans
- e) Cluster of a String and a Numeric



| ni.com/training

42

Summary—Quiz Answer

3. Which of the following are valid data types for queues?

- a) String**
- b) Numeric**
- c) Enum**
- d) Array of Booleans**
- e) Cluster of a String and a Numeric**



| ni.com/training

43

Summary—Quiz

4. The Event structure handles only one event each time it executes.
- a) True
 - b) False



| ni.com/training

44

Summary—Quiz Answer

4. The Event structure handles only one event each time it executes.
- a) True**
 - b) False



| ni.com/training

45

Lesson 2

Implementing Design Patterns

TOPICS

- A.Design Patterns
- B.Simple Design Patterns
- C.Multiple Loop Design Patterns
- D.Error Handlers
- E.Generating Error Codes and Messages
- F.Timing a Design Pattern
- G.Functional Global Variable Design Pattern



ni.com/training

46

A. Design Patterns



ni.com/training

47

Why Use Design Patterns?

- They have proven themselves useful for developing software.
- You don't have to start a program from scratch.
- They make it easier for others to read and modify your code.



Design Patterns – Code implementations and techniques that are solutions to specific problems in software design

Design patterns typically evolve through the efforts of many developers and are fine-tuned for simplicity, maintainability, and readability.



| ni.com/training

48

B. Simple Design Patterns

Simple VI

General VI

State Machine

Event-Based State Machine

Simple State Machine




| ni.com/training

49



Simple VI Pattern

- Single VI that takes a measurement, performs calculations, and either displays the results or records them to disk.
- Usually does not require a specific start or stop action from the user.




50

General VI Pattern

This pattern has three phases:


- Start-up
- Main application
- Shut-dov



ni.com/training

51

General VI Framework




ni.com/training

52

State Machine Pattern


Usually has a start-up and shut-down state, but also contains other states.



ni.com/training

53

State Machine Framework



ni.com/training

54

Event-Based State Machine

- Combines event programming with a State Machine design.
- Includes a "Wait on Event" case to processes user-interface events.




ni.com/training

55

Simple State Machine Template

Use the **Create Project** dialog box to expedite implementing an event-based state machine application.



ni.com/training

56

C. Multiple Loop Design Patterns


Producer/Consumer (Events)



ni.com/training

59

Producer/Consumer Design Patterns



ni.com/training

60

D. Error Handlers

Examples

Producer/Consumer Error Handler




ni.com/training

62

Examples of Error Handlers

- Simple Error Handler VI
- State machine error handler
- I/O error handler



ni.com/training

63

Producer/Consumer Error Handler


- Both producer and consumer loops gracefully stop when an error occurs.
 - Producer loops pass error information to consumer loops.
 - Consumer loops send stop information to producer loops.
- Transition to a Shutdown case in the consumer loop to execute shutdown code prior to stopping the VI.
- Report error information to the user.



ni.com/training

64

Error Handler VI



ni.com/training

65

E. Generating Error Codes and Messages

Error Reporting Options

Error Ring



ni.com/training

67



Error Reporting Options

• Use existing error reporting mechanisms to report error conditions detected with your code. Errors to include:

- Invalid inputs to subVIs
- File and resource errors
- LabVIEW-generated messages

• Options for error reporting:

- Pre-defined errors
- User-defined errors
- Overriding LabVIEW-generated messages



ni.com/training

68

F. Timing a Design Pattern

Execution Timing

Software Control Timing

Synchronized Timeout

Get Data/Time in Seconds



ni.com/training

72

Timing a Design Pattern

Execution Timing

- Provides the design pattern with a function that specifically allows the processor time to complete other tasks.


Software Control Timing

- Times a real-world operation to perform within a set time period.
- Controls the frequency at which a loop executes.

[ni.com/training](https://www.ibm.com/training)

73


Timing a Design Pattern

[ni.com/training](https://www.ibm.com/training)

74

Execution Timing

Design patterns where the timing is based on the occurrence of events do not need execution timing.



75

Software Control Timing

- Software control timing must allow the design pattern to run continuously without stopping.
- The following functions are better for execution timing rather than software control timing.

-Wait (ms)




-Wait Until Next ms Multiple



76

Synchronized Timeout


Even if no other events occur or the queue is empty, both the producer and consumer loops continue to execute at regular intervals.



77

Get Date/Time in Seconds

Use the Get Date/Time In Seconds function for working with elapsed times.



78

G. Functional Global Variable Design Pattern

Motivation

Problem: Read-Modify-Race Conditions

Solution: Use FGV SubVI

Function Global Variable (FGV)

Definition

Basic Framework

Other Use Cases for FGVs




ni.com/training

81

Problem: Read-Modify-Write Race Conditions


- The block diagram below has a “read-modify-write” race condition.
- Simultaneous transactions can generate incorrect data.



82

Solution: Use a SubVI Based on the Functional Global Variable (FGV) Design Pattern

Wrap the shared global resource and all functions that work on the resource into a FGV subVI.




83

The Up Down Counter VI Example



Up Down Counter VI is an FGV VI that:

- Eliminates the read-modify-write race condition for simultaneous transactions.
- Encapsulates methods for resetting, incrementing, and decrementing.





84

Up Down Counter VI

•**Action Selection**—Uses an Enum to select the counter operation:

- Reset
- Increment
- Decrement
- Get Current Value


•**Running Total**—Outputs the current value of the Up Down Counter VI.



ni.com/training

85


Up Down Counter VI Framework



ni.com/training

86

Functional Global Variable



Functional Global Variable – A non-reentrant VI that uses uninitialized shift registers to hold global data. The VI often allows for actions to be performed on the data.

- Functional**

- Can perform calculations or manipulate data.

- Global**

- Make data globally available within an Application Instance.

- Variable**

- Stores data.

- Can be written to or read from like a variable.




ni.com/training

88

The FGV Basic Framework

The general form of a functional global variable VI includes:


- An uninitialized shift register with a single iteration While Loop.
- A Case structure.



ni.com/training

89

Other Use Cases for FGVs



ni.com/training

90

Summary—Quiz

1. Which of the following are reasons for using a multiple loop design pattern?
 - a) Execute multiple tasks concurrently
 - b) Execute different states in a state machine
 - c) Execute tasks at different rates
 - d) Execute start up code, main loop, and shutdown code



ni.com/training

93

Summary—Quiz Answer

1. Which of the following are reasons for using a multiple loop design pattern?

- a) Execute multiple tasks concurrently**
- b) Execute different states in a state machine
- c) Execute tasks at different rates**
- d) Execute start up code, main loop, and shutdown code



| ni.com/training

94

Summary—Quiz

2. Which of the following are examples of error handling code?

- a) Displays a dialog box used to correct a broken VI
- b) Generates a user-defined error code
- c) Displays a dialog box when an error occurs
- d) Transitions a state machine to a shutdown state when an error occurs



| ni.com/training

95

Summary—Quiz Answer

2. Which of the following are examples of error handling code?

- a) Displays a dialog box used to correct a broken VI
- b) Generates a user-defined error code
- c) Displays a dialog box when an error occurs**
- d) Transitions a state machine to a shutdown state when an error occurs**



| ni.com/training

96

Summary—Quiz

3. What is the default timeout value of an Event structure?

- a) 0
- b) 100 ms
- c) Never time out
- d) The input value of the Wait (ms) function that exists in the same loop as the Event structure



| ni.com/training

97

Summary—Quiz Answer

3.What is the default timeout value of an Event structure?

a)0

b)100 ms

c)Never time out

d)The input value of the Wait (ms) function that exists in same loop as the Event structure



| ni.com/training

98

Lesson 3 Controlling the User Interface

TOPICS

A.VI Server Architecture

B.Property Nodes

C.Invoke Nodes

D.Control References



| ni.com/training

99

A. VI Server Architecture

VI Server Architecture

Properties and Methods

VI Class Hierarchy



| ni.com/training

100

VI Server Architecture

The VI Server provides programmatic access to LabVIEW.

Use the VI Server to:

- Programmatically control front panel objects & VIs
- Dynamically load and call VIs
- Run VIs on a computer or remotely across a network
- Programmatically access to the LabVIEW environment and editor (Scripting)



| ni.com/training

101

Properties and Methods



Properties — Single-valued attributes of the object:
read/write, read only, write only

Properties include color, position, size, visibility, label text, and label font.



Methods — Functions that operate on the object


Methods include reinitializing values to default and exporting graph images.



ni.com/training

102


VI Server—Class Hierarchy



ni.com/training

103

VI Server—Class Hierarchy



ni.com/training

104

B. Property Nodes

Definition

Creating Property Nodes

Execution Order



ni.com/training

105

Property Nodes

•Property Nodes read and write the properties of an object. Property Nodes can:



- Change the color of a chart plot.
- Disable and enable controls.
- Get the location of a control or indicator.

•Property Nodes allow you to make these modifications programmatically.

•Use Context Help to get information about properties.

•There are two types of Property Nodes.

- Implicitly linked
- Explicitly linked



ni.com/training

106

Creating Property Nodes


Create a Property Node for a front panel object.

DEMONSTRATION

107

Execution Order

- Property Nodes can have multiple properties.
- Properties execute from top to bottom.



ni.com/training

108

C. Invoke Nodes

Definition

Control Methods

VI Methods





ni.com/training

111

Invoke Nodes

- Invoke Nodes call methods or actions on objects.
 - Get VI Version
 - Print VI panel
 - Reinitialize All to Default
- Invoke Nodes perform actions on referenced items such as VIs and controls.
- Most methods have parameters.
- Use Context Help to get information on methods.
- There are two types of control Invoke Nodes.
 - Implicitly linked
 - Explicitly linked




ni.com/training

112

Control Methods

- To create an implicitly linked Invoke Node:
 1. Right-click the control terminal on the block diagram and select **Create»Invoke Node**.
 2. Select a method from the submenu.
- Examples of control methods:
 - Reinitialize to Default
 - Export Image




ni.com/training

113

VI Methods

- Use a VI Server Reference to associate an Invoke Node with the current VI.
- To create a VI method:
 1. Place an Invoke Node on the block diagram.
 2. Right-click and select **Select Class** to choose a class.
 3. Right-click again and select **Select Method** to choose a method.



ni.com/training

114

D. Control References

Implicitly and Explicitly Linked Property Nodes

Create SubVIs

Use the Create SubVI Tool

Create Manually

Select a VI Server Class




ni.com/training


117

Control References

Implicitly Linked Property Node



Explicitly Linked Property Node




ni.com/training

118

Control References

- A control reference is a reference to a front panel object.
- Wire control references to generic Property Nodes.
- Pass control references to subVIs.





ni.com/training

119

Creating a SubVI

To create explicitly linked Property Nodes in a subVI:



1. Create your VI.
2. Select the portion of the block diagram that will be in the subVI.
3. Select **Edit»Create SubVI**.
LabVIEW automatically creates the control references needed for the subVI.
4. Customize and save the subVI.



ni.com/training

120


Creating a SubVI



ni.com/training

121

Create Control References Manually



122

Selecting the VI Server Class

- After you place a Control Refnum on the front panel of a subVI, specify the VI Server class of the control.
 - Right-click and select **VI Server Class** from the shortcut menu.
 - Alternatively, drag a control into a control refnum to specify the type.
- The class specifies the type of control references that the subVI accepts.



123

Summary—Quiz

1. For each of the following items, determine whether they operate on a VI class or a Control class.

- a. Format and Precision
- b. Visible
- c. Reinitialize to Default Value
- d. Show Tool Bar



| ni.com/training

129

Summary—Quiz Answer

1. For each of the following items, determine whether they operate on a VI class or a Control class.

- a. Format and Precision: **Control**
- b. Visible: **Control**
- c. Reinitialize to Default Value: **Control**
- d. Show Tool Bar: **VI**




| ni.com/training

130

Summary—Quiz

2. You have a Numeric control refnum in a subVI. Which control references could you wire to the control refnum terminal of the subVI?

- a. Control reference of a Knob
- b. Control reference of a Numeric Array
- c. Control reference of a Thermometer indicator
- d. Control reference of an LED




ni.com/training

131

Summary—Quiz Answer

2. You have a Numeric control refnum in a subVI. Which control references could you wire to the control refnum terminal of the subVI?

- a. Control reference of a Knob**
- b. Control reference of a Numeric Array
- c. Control reference of a Thermometer indicator**
- d. Control reference of an LED



ni.com/training

132

Lesson 4

File I/O Techniques

TOPICS

- A. Compare File Formats
- B. Create File and Folder Paths
- C. Write and Read Binary Files
- D. Work with Multichannel Text Files with Headers
- E. Access TDMS Files in LabVIEW and Excel



ni.com/training

133

A. Compare File Formats




ni.com/training

134

Compare File Formats

At their lowest level, all files written to your computer's hard drive are a series of bits.



ni.com/training

135

Compare File Formats

Numeric Precision	Good	Best	Best
Share data	Best (Any program easily)	Better (NI Programs easily; Excel)	Good (only with detailed format information)
Efficiency	Good	Best	Best
Ideal Use	Share data with other programs when file space and numeric precision are not important.	Store measurement data and related metadata. High-speed streaming without loss of precision.	Store numeric data compactly with ability to random access .



ni.com/training

136

B. Create File and Folder Paths

Methods of Creating

Creating Relative Paths and Folders

Dynamically Creating Filenames



| ni.com/training

137

Methods of Creating File and Folder Paths

- Hard-coded paths
 - Useful for quick prototypes.
 - Not recommended for applications.
- File Dialog
 - Allow user to specify the path to a file or directory.
 - Customize dialog options to limit options (*.txt).
- Programmatic creation
 - Create consistent filenames and extensions.
 - Example: testdata_001.txt, testdata_002.txt, etc.
 - Specify a consistent location.




| ni.com/training

138

Creating Relative Paths


- Relative paths set paths relative to the application or system directory.
- Use the Application Directory VI to get project directory paths.
- Use the Get System Directory VI to get system directory paths.
- Paths differ based on the operating system and user.



139


Creating Folders

Check if folder already exists before creating. Otherwise, attempting to create an existing folder will result in an error.



140

Dynamically Creating Filenames



ni.com/training

141

C. Write and Read Binary Files

Using Binary File Functions
Demonstration




ni.com/training

144

Using Binary File Functions

Use Binary File functions to read and write binary files.

- You can create custom file types.
- You must know the required file formats for your system.



145

Demonstrate Writing a Bitmap File

Use Binary File I/O to write a bitmap format file.

DEMONSTRATION

146

D. Work with Multichannel Text Files with Headers

Review of Text Files

Add Headers to the File

Write Multichannel Data

Read Data and Extract Information



| ni.com/training

147

Review of Text Files

- Use ASCII characters to store information in text files.
- Each character (number, letter, punctuation) takes 1 byte.
- Many applications can open text files, including Excel.
- Files are typically larger and slower to write/read than binary.
- You cannot randomly access text files.



| ni.com/training

148

Creating Text Files with Headers

No Header Data

24.45		34.54	
23.41		35.32	
22.97		35.98	
21.56		36.76	



Tab




CR/LF

Header Data

Operator Name	David		
UUT S/N	A1234		
Test Name	Pressure		
Channel Name	Temperature	Pressure	
Units	Kelvin	PSI	
Max. Value	24.45	36.76	
	24.45	34.54	
	23.41	35.32	
	22.97	35.98	
	21.56	36.76	

149


Creating Text Files with Headers— Hard Coding



Hard coding header information becomes difficult to maintain.

150

Creating Text Files with Headers—SubVIs



- Using subVIs to build strings is a scalable solution.
- With subVIs you can add and remove header data easily.
- The block diagram is more readable.




ni.com/training

151

Writing Multiple Channels

- LabVIEW stores multi-dimensional arrays in row-major order.
 - Rows are identified by the first index of a 2-D array.
 - Columns are identified by the second index.
- Transpose data before writing to file to view channel data in column format.




ni.com/training

152

Reading Channel Data

- Use Array Subset to skip over header information.
- Transpose data after reading to convert data back to row-major order.
- Use Index Array to extract one column or row of data.



153

Read Data and Extract Information


- After reading a spreadsheet file, how do you extract information?
- How do you find a property value?
- How do you read the UUT Serial Number value?
- How do you extract channel data for column 1?

Operator Name	Fred	
UUT Serial Number	A001	
Test Name	Stress Response	
Channel Data		
0	0.000000	0.000000
1	0.049068	0.309017
2	0.098017	0.587785
3	0.146730	0.809017
4	0.195090	0.951057
5	0.242980	1.000000

156

Given a Property Name, Find the Value

- Use Index Array to put column 0 in a 1-D array.
- Use Search 1D Array to find Property Name:
 - Return row index if found.
 - Return -1 if not found.
- Use Index Array to get element at row index and column 1.
- What should you do if Property Name is not found?




ni.com/training

157

Extracting a Data Channel

- Identify row index of first data value.
- Use Array Subset to extract Data Channel.
 - Output is one channel of data in a 2-D string array.
- Convert to 1-D numeric array.



ni.com/training

158

E. Access TDMS Files in LabVIEW and Excel

TDMS File Format

Data Hierarchy and Properties

TDMS Functions

File Viewer



| ni.com/training

159

TDMS File Format

The Technical Data Management Streaming (TDMS) file format contains two types of data:

- Meta data – Names and properties
- Raw data – Measurement data in binary format

Use TDMS files for the following purposes:

- To store test or measurement data.
- To create structures for grouping your data.
- To store headers/properties about your data.
- To read and write data at high speeds.


A variety of applications, including Microsoft Excel, can access TDMS files.



| ni.com/training

160


TDMS Files—Data Hierarchy and Properties



ni.com/training

161

TDMS Functions




ni.com/training

162

TDMS Files—File Viewer

Open TDMS files and present the file data in the TDMS File Viewer dialog box.



163

Exercise 4-3 Write and Read TDMS Files

Learn how to read data from a TDMS file.


GOAL

164

Summary—Quiz

1. Consider the following code. The resulting Log File Path contains a text file path in which folder?

- a) Same folder as VI that executed the code.
- b) Same folder as the LabVIEW project.
- c) Current user's AppData directory.
- d) Same folder as the Application Directory VI.



ni.com/training

166

Summary—Quiz Answer

1. Consider the following code. The resulting Log File Path contains a text file path in which folder?

- a) Same folder as VI that executed the code
- b) Same folder as the LabVIEW project**
- c) Current user's AppData directory
- d) Same folder as the Application Directory VI




ni.com/training

167

Summary—Quiz

2. What index value is returned from Search 1D Array function if Property Name is not found in the input array?

- a) NaN (Not a Number)
- b) 0
- c) -1
- d) Negative Infinity




ni.com/training

168

Summary—Quiz Answer

2. What index value is returned from Search 1D Array function if Property Name is not found in the input array?

- a) NaN (Not a Number)
- b) 0
- c) -1
- d) Negative Infinity



ni.com/training

169

Summary—Quiz

3. You need to store data that other engineers will later analyze with Microsoft Excel. Which file storage format(s) should you use?

- a) Tab-delimited ASCII
- b) Custom binary format
- c) TDMS



| ni.com/training

170

Summary—Quiz Answer

3. You need to store data that other engineers will later analyze with Microsoft Excel. Which file storage format(s) should you use?

- a) Tab-delimited ASCII**
- b) Custom binary format
- c) TDMS**



| ni.com/training

171

Summary—Quiz

4. TDMS files store properties at which of the following levels?

- a) File
- b) Channel Group
- c) Channel
- d) Value



| ni.com/training

172

Summary—Quiz Answer

4. TDMS File store properties at which of the following levels?

- a) File**
- b) Channel Group**
- c) Channel**
- d) Value



| ni.com/training

173

Lesson 5

Improving an Existing VI

TOPICS

- A. Refactoring Inherited Code
- B. Typical Refactoring Issues



ni.com/training

174

A. Refactoring Inherited Code

Definition

When To Refactor

The Refactoring Process



ni.com/training


175

Refactoring Inherited Code

Inherited VIs may be poorly designed, making it difficult to add features later in the life of the VI.

Refactoring:

- Is the process of redesigning software to make it more readable and maintainable so that the cost of change does not increase over time.
- Changes the *internal* structure of a VI to make it more readable and maintainable, without changing its observable behavior.



ni.com/training

176

When to Refactor

- When you are adding a feature to a VI or debugging it.
- Good candidates for complete rewrites:
 - VIs that do not function.
 - VIs that satisfy only a small portion of your needs.



ni.com/training

177

Refactoring Process

When you refactor to improve the block diagram, make small cosmetic changes before tackling larger issues.



| ni.com/training

178

B. Typical Refactoring Issues

Disorganization
Poor Naming
Overly Complicated
Duplicated Logic
Breaks Dataflow
Outdated Practices



| ni.com/training

179

Disorganization

The block diagram is too disorganized, too big, or includes too many nested structures

Solution:

- Move objects within the block diagram.
- Create subVIs to make it smaller and more organized.
- Place comments to improve readability.




ni.com/training

180

Poor Naming

The block diagram uses incorrect object names and poor icons.



ni.com/training

181


182

Duplicate Logic

The block diagram uses duplicate logic.

Solution:

Refactor the VI by creating a subVI for the duplicated logic.




183

Breaks Dataflow

The block diagram does not use dataflow programming.

Solution:

- Replace Sequence structures with state machines.
- Delete local variables and wire directly to controls or indicators.




184

Outdated Practices

The VI was created in an earlier version of LabVIEW and has outdated practices.

Solutions:

- Replace polling-based design with event-based design.
- Use new features that simplify code.



185

Summary—Refactoring Checklist

Use the following refactoring checklist to help determine if you should refactor a VI:

- ☐ VI is too disorganized, too big, or includes too many nested structures.
- ☐ VI uses incorrect object names and poor icons.
- ☐ VI uses unnecessary logic or has complicated algorithms.
- ☐ VI has duplicated logic.
- ☐ VI does not use dataflow programming.
- ☐ VI uses outdated development practices.



ni.com/training

188

Lesson 6 Creating and Distributing Applications

TOPICS

- A. Preparing the Files
- B. Build Specifications
- C. Create and Debug an Application
- D. Create an Installer



ni.com/training

189

A. Preparing the Files

Building Applications Checklist

VI Properties

Paths and System Paths

Quit Application



| ni.com/training

190

Preparing the Files

To create a professional, stand-alone application for your VIs, you must first prepare your files:

- Recompile and save changes to VIs.
- Verify desired VI Properties settings.
- Ensure paths generate correctly.
- Conditionally call the Quit LabVIEW function.



| ni.com/training

191

Building Applications Checklist

Refer to the *Building Applications Checklist* topic of the *LabVIEW Help* for more information about:

- Preparing files.
- Configuring build specifications.
- Distributing builds.





ni.com/training

192

Preparing the Files – VI Properties

To set VI Properties settings:

- Manually edit the **VI Properties** dialog box.
- Programmatically change VI Properties using VI Server.




ni.com/training

193

Preparing the Files – Paths

Set paths relative to Application Directory VI path.

From a stand-alone application	The path to the folder containing the application executable
From a LabVIEW project (.lvproj) in the LabVIEW Development Environment	The path to the project folder



ni.com/training

194

Preparing the Files – System Paths

- Use the Get System Directory VI to get system directory paths.


- Paths differ based on the operating system and user.

- Windows XP User Documents:

-C:\Documents and Settings\<user>\My Documents

- Windows 7 User Documents:

-C:\Users\<user>\Documents




ni.com/training

195

Preparing the Files – Quit Application

Programmatically quit your application using the Quit LabVIEW Function.

To conditionally call the Quit LabVIEW function, use the App.Kind property.



ni.com/training

196

Exercise 6-1A Preparing Files for Distribution

- Review the Building Applications Checklist.
- Prepare VIs to build a standalone application.

GOAL

197

Exercise 6-1A

Preparing Files for Distribution

What are some of the changes you observed when you ran your VI after changing the Window Appearance properties?

DISCUSSION

198

B. Build Specifications

199

Build Specifications

A build specification contains all the settings for the build, such as files to include, directories to create, and settings for VIs.




ni.com/training

200

Why Use Build Specifications?

Use build specifications to build the following:

- Stand-alone applications
- Installers
- Source distributions
- Zip files
- Shared libraries
- Packed Project libraries
- .NET Interop Assemblies
- Web services



201

C. Create and Debug an Application

System Requirements for Applications

Configure the EXE Build Specification

Build and Run the EXE

Debug the EXE



| ni.com/training

202

System Requirements for Applications

- Applications that you create with Build Specifications generally have the same system requirements as the LabVIEW development system used to create the VI or application.
- Memory requirements vary depending on the size of the application created.



| ni.com/training

203

Configure the EXE Build Specification

In the **Application Properties** dialog box:


- Specify name of executable.
- Specify known destination for generated executable files
- Identify a startup VI and include any dynamically linked files.
- Enable debugging, if desired.
- Preview the build.



ni.com/training


204

EXE Properties – Information




205

EXE Properties – Source Files



206

EXE Properties – Advanced



207

208



| ni.com/training


209



| ni.com/training

210


Run EXE and Verify Execution



211

212

Debug the EXE from LabVIEW



213

Why Executables Might Behave Differently

- File paths may change, which can lead to errors.
- Drivers or support files are missing.
- System resources like memory or CPU speed may differ, which can lead to timing changes.
- The application INI file differs from the LabVIEW INI file.
- Not all features are supported in the Run-Time Engine.



| ni.com/training

214

D. Create an Installer

Why Create an Installer?

Configure the Installer Build Specifications

Deploy the Application to Another Machine

Debug the Executable



| ni.com/training

217

Why Create an Installer?

- Executables need the LabVIEW Run-Time Engine (RTE) to execute on a target system.
- If an application relies on drivers, those drivers need to be installed on the target system.
- An installer ensures files are copied to the right places.
- Professional applications use installers.



| ni.com/training

218

Configuring Installer Build Specifications

In the **Installer Properties** dialog box:

- Include the LabVIEW Run-Time Engine.
 - From the Additional Installers page, select the **NI LabVIEW Run-Time Engine**.
- Include drivers used in the application.
 - For example, if you use DAQmx VIs in your application, then include the NI-DAQmx driver.


Refer to the *Caveats and Recommendations for Building Installers* topic of the *LabVIEW Help* for more information.



| ni.com/training


219

Installer Properties – Product Information




220

Installer Properties – Additional Installers



221

Install and Run on Destination Computer



222

Debug Executable on Destination Computer

1. Find IP address of destination computer.
 - Type “`ipconfig`” in command prompt window of the destination computer to discover IP address.
2. Launch LabVIEW on development machine.
3. Enter the IP address or machine name in **Debug Application or Shared Library** dialog box.
4. Click the **Connect** button in the **Debug Application or Shared Library** dialog box.
 - Control is transferred to development machine.
5. Debug the executable using LabVIEW debugging tools.

223

Summary

- The Application Builder enables you to create stand-alone applications and installers.
- To create a professional, stand-alone application you must understand:
 - The architecture of your application.
 - The programming issues particular to your application.
 - The application building process.
 - The installer building process.
- Test your application and installer regularly as you develop.



| ni.com/training

226

Thank you!



| ni.com/training

235